


QTA Tracer System® - Overview



1 Tracer settings are configured in **QTA Access Point Configurator**.

2 Tracer is initiated (linked to specific blood bag) and scanned for status in **QTA Access Point**. Information uploaded to **QTA Web Portal**.

3 Check Status at any time by turning blood bag with tracer in your hand. Green light - blood ok. Red light - expiration date passed.

4 Blood Status, Temperature logs, location etc - easily accessible on the **QTA Web Portal**.

5 **QTA Check-in Nodes** (iPad with free application) strategically placed throughout hospitals (and accompanying transport) register tracers passing nearby and uploads information to the QTA Web portal

6 Scan tracer with **QTA Reader** (Smartphone with free application) for blood status and expiration date.

